

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

SHANGHAI REYOUNGEL MEDICAL TECHNOLOGY COMPANY LIMITED

Add: Floor 15, Building 2#, No.1809, Qixin Road, Minhang District, 201100, Shanghai, China
Tel: 0086 21 51698555 E-mail: info@reyoungel.com Web: www.reyoungel.com

See What
We Can Do For Your
Beauty

COMPANY INTRODUCTION

Bioha Laboratories was established in 2003, who created the Reyoungel brand for dermal filler. Shanghai Reyoungel Medical

Technology Company Limited is the branch for Reyoungel brand, which has more than 15 years history in manufacturing and exporting Reyoungel® dermal filler. Our facility is ISO 13485

This is a watermark for the trial version, register to get the full one!

We are a professional company focusing on the development of dermal and beauty products. Our products range includes

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

During the aging process, the human face loses fat and volume while the skin loses collagen and elasticity. Reyoungel® products are developed to smooth outlines and wrinkles, shape facial contours and add volume to lips. It is very safe and effective to treat the entire face creating a more beautiful and youthful look.

Treatment Areas

Smooth forehead furrows

Fill in frown lines

Remove crow's feet

Remove bags and dark circles under eyes

Facial sculpt, and smooth cheek

Smooth marionette lines

Contour chin and jaw

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

Products Range



Fine Lines 1ml/2ml

Ingredient: Cross-linked hyaluronic acid gel

HA Concentration: 20mg/ml

Particle Size: 0.10-0.15mm

Duration: 6-9 months

Function:

Fill superficial lines and wrinkles, crow's feet, perioral lines, oral commissures, etc.



Derm Plus 1.5ml/10ml

Ingredient: Cross-linked hyaluronic acid gel

HA Concentration: 20mg/ml

Particle size: 0.5-1.25mm

Duration: 9-18 months

Function:

Derm plus is designed to fill deep wrinkles, recover and form volume of cheeks, cheekbones and chin, body contouring.

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now



Derm Deep 1ml/2ml/10ml

Ingredient: Cross-linked hyaluronic acid gel

HA Concentration: 20mg/ml

Particle Size: 0.28-0.5 mm

Duration: 6-12 months

Function:

Fill deep wrinkles, like nasolabial folds. Shape facial contour, lips and cheek volume expansion.

Products Info

	Fine Lines	Derm	Derm Deep	Derm Plus	Sub Skin
Concentration	20mg/ml	20mg/ml	20mg/ml	20mg/ml	20mg/ml
Particles Size (mm)	0.10-0.15	0.15-0.28	0.28-0.5	0.5-1.25	1.25-2.0
Recommendation	Thin superficial lines, such as worry lines, forehead lines, and crow's feet.	Moderate wrinkles, such as glabellar, oral commissures, and lip corners.	Deep facial wrinkles and folds, like the nasolabial folds, shaping cheeks, chin, and lips augmentation.	Fill deep wrinkles, recover and form volume of cheeks, chin, and lips contouring.	Form volume, correct facial contours, expand volume of cheeks, chin, restore and contouring.
Inject					
BDDF	< 0.2 µg/g	< 0.2 µg/g	< 0.2 µg/g	< 0.2 µg/g	< 0.2 µg/g
Packaging (pcs)	2ml	2ml	2ml/10ml	10ml	10ml

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:


- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

How it works?



Hyaluronic acid is the key component of the skin to hold water. As we age, the amount of hyaluronic acid in the skin decreases, leading to wrinkles and fine lines. Hyaluronic acid is a naturally occurring polysaccharide in the skin; it's considered to have a greater capacity to hold moisture. In fact, one gram of hyaluronic acid is able to hold up to six liters of water. The skin is generally comprised of 3 major layers, the epidermis, dermis and subcutis. Hyaluronic acid exists in all layers of the skin. It adds moisture to the skin, and help it retain that moisture. When hyaluronic acid is injected into the dermal layer, they can have the effect of temporarily plumping wrinkles and fine lines.


Remove Watermark Now



- Fine · Upper part of dermis
- Derm · Middle part of dermis
- Deep · Deep layer of dermis
- Subskin · Subcutaneous tissue

Study Datas


Extrusion Force



This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.



Injection Reaction




Remove Watermark Now

	Standard	Result
Appearance	Colorless, transparent, foreign without impurities	Meet specification
Content	16.0-24.0mg	20.3mg
pH	6.5-7.5	7.1
BDDE* Residues	≤0.2μg/ml	< 0.2μg/ml
Heavy Metal	≤10μg/ml	<10μg/ml
Ethanol	≤0.15%	0.011%

Conclusion: The product conforms to the standard of Enterprise Standard.

*What is BDDE?
It is a chemical cross-linking agent that produces Hyaluronic Acid into a gel form.

Testimony



· Frown lines filled with Reyoungel® Fine Lines

Safe & Natural

Reyoungel® is from non-animal origin and high-purity, it can be well biodegraded in the skin. There is no side effect and less invasive, filling wrinkles and smoothing your skin without inhibiting facial expression.

Fast & Effective

Patients need no hospitalization and surgery, the whole process takes only 30 minutes, and there are obvious results after the first treatment.

This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now



· Lip volumized with Reyoungel® Derm

93%

Doctors Satisfied

Doctors consider Reyoungel® is easy to inject, and has its complete range products designed to match different level to the tissue structures.

87%

Patients Satisfied

Reyoungel® can effectively enhance the natural looks by restoring skin's volume, bring patients a nature, safe and effective beauty.



This is a watermark for the trial version, register to get the full one!

Benefits for registered users:

- 1.No watermark on the output documents.
- 2.Can operate scanned PDF files via OCR.
- 3.No page quantity limitations for converted PDF files.

Remove Watermark Now

